

Low and Ultra-Low Positive COVID-19 Controls

Thermo Scientific AcroMetrix Coronavirus 2019 (COVID-19) RNA Control (RUO)


The Thermo Scientific™ AcroMetrix™ Coronavirus 2019 (COVID-19) RNA Control was designed and developed as a multiplex synthetic RNA transcript controls used to validate and monitor the performance of molecular tests for detecting COVID-19. This is a positive control kit containing two concentrations: a low positive and ultra-low positive concentration to avoid false-negative results. Both the controls are quantitated using Bio-Rad™ Droplet Digital™ PCR (ddPCR™) technology.

Product Features


- Multiplex control containing highly unique N, S, E and ORF1ab genome regions of SARS-CoV-2
- Product stable up to 3 freeze/thaws
- Open vial stability at 2-8°C for 7 days
- Product storage -20°C or below

AcroMetrix Coronavirus 2019 (COVID-19) RNA Control sequence design is compatible with, but not limited to, the following assays:

- TaqPath™ COVID-19 Combo kit (A47814)
- USA CDC assay
- China CDC assay
- Japan National Institute of Infectious Diseases assay
- Thailand National Institute of Health assay
- Charité – Universitätsmedizin Berlin assay



For Research Use only



Catalog Number	Control Name	Copies/μL	Quantity
954519	AcroMetrix COVID-19 Low Positive Control	500	1 x 20 μL
	AcroMetrix COVID-19 Ultra-Low Positive Control	100	1 x 20 μL

Customer Support

Phone: 1-800-232-3342 - Press 1 for Customer Service, then 3 for QC

Fax: 1-800-829-8115

e-mail: mgc-masorders@thermofisher.com

Sales Support

Phone: 1-800-232-3342

e-mail: sales.diagnostics.fmt@thermofisher.com

Technical Service and Support

Phone: 1-800-232-3342- option 2, then 2 for QC

Clinical Diagnostics

USA

46500 Kato Road
Fremont, CA 94538 USA
Tel. 800-232-3342 / + 1 510-979-5000
sales.diagnostics.fmt@thermofisher.com

Asia Pacific

Tel: +65 6499 9999
cdd.asia.info2@thermofisher.com

Australia

Tel. + 61 1800 333 110
auiinfo@thermofisher.com

Canada

Tel. 905-286-4290
CDD.Canada@thermofisher.com

China

Tel. +86 400 650 5118
sales.china@thermofisher.com

France, Belgium (Brussels / Wallonia) & Luxembourg

Tel. + 33 (0) 1 40 86 65 00
cdd.fr.info@thermofisher.com

Germany

Tel. +49 (0) 3302 883 608
cdd.de.info@thermofisher.com

Hong Kong

Tel. + 852 3107 7600
sales.hk@thermofisher.com

Italy

distributor.cdd@thermofisher.com

Japan

Tel. + 81 (0)120-147-075
JPYOK-CDD.QC@thermofisher.com

Netherlands & Belgium (Flemish)

Tel. +49 3302 883 660
Info.benelux.cdd@thermofisher.com

New Zealand

Tel. + 61 0800933 966
nzinfo@thermofisher.co.nz

Nordic

Tel. + 46 8 594 691 30
info.nordic.cdd@thermofisher.com

Russia

Tel. + 7 495 739 76 41
info.btd.moscow@thermofisher.com

South Africa

Tel. +27 11 792 6790
support-za.idd@thermofisher.com

Spain & Portugal (Iberia)

Tel. + 34 93589 8338
CDD.Info.ES@thermofisher.com

Switzerland & Austria

Tel. + 41 26 663 86 70
cdd.ch.info@thermofisher.com

United Kingdom & Ireland

Tel. + 44 1442 868 940
cdd.UK.info@thermofisher.com

For Countries Not Listed

distributor.cdd@thermofisher.com

Learn more at: [thermofisher.com/covid-19-control](https://www.thermofisher.com/covid-19-control)

Availability of product in each country depends on local regulatory marketing authorization status.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. Bio-Rad, Droplet Digital, and ddPCR are trademarks of Bio-Rad. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

FR-MTL-10027916-0-AcroMetrix-COVID-19-Flyer-EN

ThermoFisher
SCIENTIFIC